

# XPE Accumulator Patient Table



## Digital Radiography of Maximum Flexibility

Off-grid, mobile patient table with motorized elevating, floating table top. The easy usage allows an optimal patient positioning. The table is stand-sided open and makes a positioning right over the digital detector possible.



Data sheet Patient Table

## Technical Data

### Table Top

Table top length	200 cm
Width	75,5 cm
Height	58 – 84 cm
Transverse movement	20 cm
Longitudinal movement	92 cm
Table top thickness	2,5 cm
Table top material	Carbon fibre
Alu equivalent	0,7 mm AL
Weight	98 kg
Max. patient weight	150 kg
Accumulator charge station	Power connection 100 – 240 V, 50 / 60 Hz, 1,2 A incl. 2 pcs. accumulator

### Standardization

Mechanical safety	Acc. to IEC/EN 60601-1, UL 2601-1
X-ray	Acc. to IEC/EN 60601-1-3

### Accessory

Produced by Provotec  
Handgrip



Gesellschaft für  
medizinische Datenverarbeitung  
und Integration mbH  
Schuhstr. 30  
D-91052 Erlangen

**Tel:** + 49 (0) 91 31 - 6 90 87 - 40  
**Fax:** + 49 (0) 91 31 - 6 90 87 - 50  
eMail: [info@medigration.de](mailto:info@medigration.de)  
web: [www.medigration.de](http://www.medigration.de)



For more information please visit  
[www.medigration.de](http://www.medigration.de)  
If you have any further questions  
please feel free to contact us!