





Patient-CD Importer

Fast and simple import of patient CDs

- Intuitively operated software for the import of patient CDs
- Quick overview of the image data due to ergonomic interface
- Simple option for changing patient name and ID
- Convenient import of PDF files and findings into the PACS archive
- Optional assistance for connection to a DICOM work list



DIGITAL IMAGES ARE OUR WORLD

The medigration Patient-CD Importer is an intuitively operated piece of software that allows you to import patient CDs in just a few steps, with no problem or risk. The ergonomic user interface allows a swift overview of the respective image data. It can be used to specifically select individual images/studies or to import all images with a click. As well as patient CDs, PDF files can also be swiftly and conveniently saved into the PACS archive using the CD importer.

Function

The medigration Patient-CD Importer not only allows you to save patient data into your archive safely and without problems, but also offers you the option of asking whether a patient in your work list has been registered for an examination. For the work list, as well as for the import into the archive and searching through the CD you have at your disposal search fields for the name, ID, date of birth and sex of the patient, so that you can find the right patient and assign the images and PDF files to them.

Furthermore, in the overview section, you have the option of double-clicking on an image, loading it and displaying it in a window, zooming into it or shifting it. This makes it easier for you to pre-select the images and series before importing.

Requirements

The Patient-CD Importer can be installed on any workplace computer with a Windows operating system (2000/XP/Win 7). Regardless of whether the image data is imported straight away at front-desk registration or at a findings workplace, the user-friendly software is all-purpose and runs without any problems at all workplaces.



Ergonomic user interface for swift orientation



Gesellschaft für
medizinische Datenverarbeitung
und Integration mbH
Schuhstr. 30
D-91052 Erlangen

Tel: + 49 (0) 91 31 - 6 90 87 - 40
Fax: + 49 (0) 91 31 - 6 90 87 - 50
eMail: info@medigration.de
web: www.medigration.de



For more information please visit
www.medigration.de.
If you have any further questions
please feel free to contact us!